

## SPECIALISED CALL DEVICES

### SUMMARY

Advance Care's Specialized Nurse Call Devices have been specifically designed to assist residents with varying degrees of movement impairment. All devices come equipped with the industry standard fully moulded 6.35mm jack.



#### Head Activated Call Device

The Head Activated Nurse Call Device is a sensitive touch pad, which can be utilized in a variety of environments, and a number of different applications.

This device can be positioned to be activated by the foot, hand, elbow or head, thus ensuring ease of usage for residents with many different movement limitations.

The Head Activated Nurse Call Device comes with two secure clips and a Velcro strap. It is easily mounted onto a bed rail, or attached to the resident's pillow or clothing.

#### Breath Activator

The Breath Activator device is designed for those residents with severely limited movement.

The unit can be securely clamped to any adjacent fitting, such as a bed rail, and the flexible adjustment arm allows the device to be positioned in the most appropriate location to best suit the residents needs.

The filtered mouth piece on the unit is replaceable. The Breath Activator is UL 1069 listed and is able to be used in oxygen rich environments.

#### Geriatric Call Pendant

The Geriatric Call Pendant only requires limited pressure for activation and is suitable for a wide range of health environments. The pad is covered in highly resilient rubberized moulding.

The Pendant can be securely clipped to bedding or garments, to ensure it stays within the resident's reach.

| PRODUCT CODE | DESCRIPTION                       |
|--------------|-----------------------------------|
| AC - BA      | Breath Activated Call Device      |
| AC - BA / MP | Mouth Piece replacement for above |
| AC - GCC     | Geriatric Call Pendant            |
| AC - HA      | Head Activated Call Device        |